

# Heartsaver® Infant CPR Skills Testing Checklist



Student Name \_\_\_\_\_ Date of Test \_\_\_\_\_

Scenario: “While you are pushing a baby in a stroller at the park, you notice something is wrong with the baby. You do not have a phone nearby. You ensure that the scene is safe and take the baby out of the stroller. Demonstrate what you would do next.”

## Assessment and Activation

- Checks responsiveness     Shouts for help     Checks breathing

Once student shouts for help, instructor says, “No one is around to help.”

## Cycle 1 of CPR (30:2) \*CPR feedback devices preferred for accuracy

### Infant Compressions

- Performs high-quality compressions\*:
- Placement of 2 fingers in the center of the chest, just below the nipple line
  - 30 compressions in no less than 15 and no more than 18 seconds
  - Compresses at least one third the depth of the chest, about 1½ inches (4 cm)
  - Complete recoil after each compression

### Infant Breaths

- Gives 2 breaths with a barrier device:
- Each breath given over 1 second
  - Visible chest rise with each breath
  - Gives 2 breaths in less than 10 seconds

## Cycle 2 of CPR (repeats steps in Cycle 1) Only check box if step is successfully performed

- Gives 30 high-quality compressions     Gives 2 effective breaths

## Cycle 3 of CPR (repeats steps in Cycle 1) Only check box if step is successfully performed

- Gives 30 high-quality compressions     Gives 2 effective breaths

Instructor says, “You have just completed 5 sets of 30 compressions and 2 breaths.”

## Activates Emergency Response System (9-1-1)

- Verbalizes the need to leave to phone 9-1-1

**STOP TEST**

### Instructor Notes

- Place a ✓ in the box next to each step the student completes successfully.
- If the student does not complete all steps successfully (as indicated by at least 1 blank check box), the student must receive remediation. Make a note here of which skills require remediation (refer to Instructor Manual for information about remediation).

**Test Results**    Check **PASS** or **NR** to indicate pass or needs remediation:

**PASS**    **NR**

Instructor Initials \_\_\_\_\_ Instructor Number \_\_\_\_\_ Date \_\_\_\_\_

# Heartsaver®

## Infant CPR

### Skills Testing Critical Skills Descriptors

- 1. Assesses infant and activates emergency response system (this *must* precede starting compressions) within a maximum of 30 seconds. After determining that the scene is safe:**
  - Checks for responsiveness by tapping and shouting
  - Shouts for help
  - Checks for no breathing or no normal breathing (only gasping)
    - Scans from the head to the chest for a minimum of 5 seconds and no more than 10 seconds
- 2. Cycle 1: Performs high-quality chest compressions (initiates compressions immediately after recognition of cardiac arrest)**
  - Correct placement of fingers in center of chest
    - 2 fingers just below the nipple line
  - Compression rate of 100 to 120/min
    - Delivers 30 compressions in 15 to 18 seconds
  - Compression depth and recoil—compress at least one third the depth of the chest, about 1½ inches (4 cm)
    - Use of a commercial feedback device/manikin is highly recommended
    - Complete chest recoil after each compression
- 3. Cycle 1: Provides 2 breaths by using a barrier device**
  - Opens airway adequately
    - Uses a head tilt–chin lift maneuver
  - Delivers each breath over 1 second
  - Delivers breaths that produce visible chest rise
  - Avoids excessive ventilation
  - Resumes chest compressions in less than 10 seconds
- 4. Cycle 2: Performs same steps for compressions and breaths as in Cycle 1**
- 5. Cycle 3: Performs same steps for compressions and breaths as in Cycle 1**